ABSORBABLE SURGICAL SUTURE CATGUT PLAIN

Type Twisted multifilament with a monofilament appearance

Composition Strands of purified collagen taken from the serosal layer of selected bovines.

Origin from BSE-free classified countries.

COATING None

Tissue reaction Moderate

Absorption By phagocytosis. Complete mass absorption in approximately 63 days.

Approximately 50% of tensile strength remain after 7 days. When used in infected tissues or in tissues with increased levels of proteolytic enzymes

(stomach, cervix, vagina) it is more quickly absorbed.

Tensile strength Post implantation Approximate % original strength remaining

7 days 50 % 14 days 20 %

SUTURE COLOUR IVORY

RANGE USP 5/0 - 2/0

INDICATIONS Catgut plain sutures are intended for use in general soft tissue closing and/or

ligation; especially in general surgery, gastrointestinal surgery, gynaecology,

obstetrics, urology, ophthalmic surgery









	Round bodied needle – taper point – 1/2 circle											
NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER				
HR37	37 mm		3.5	2/0	75 cm	ivory		CP20HR37				

Reverse cutting needle – 3/8 circle												
NEEDLE			EP	USP	LENGTH	SUTURE COLOUR	OTHER INFORMATION	ORDER NUMBER				
DS16	16 mm		1.5	5/0	75 cm	ivory		CP50DS16				
			2	4/0	75 cm	ivory		CP40DS16				
DS19	19 mm		3	3/0	75 cm	ivory		CP30DS19				
DS24	24 mm		3.5	2/0	75 cm	ivory		CP20DS24				